

# Recombinant Human Brain Natriuretic Peptide/BNP Protein (His & Flag Tag)

Catalog No. PKSH032127

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

## **Description**

Synonyms Natriuretic peptides B;Gamma-brain natriuretic peptide;NPPB;BNP

SpeciesHumanExpression HostE.coli

Sequence His27-Arg102

AccessionP16860Calculated Molecular Weight11.0 kDaObserved molecular weight14 kDaTagN-His-Flag

**Bioactivity** Not validated for activity

## **Properties**

**Purity** > 95 % as determined by reducing SDS-PAGE.

**Endotoxin** < 1.0 EU per µg of the protein as determined by the LAL method.

**Storage** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to

-80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots

of reconstituted samples are stable at < -20°C for 3 months.

**Shipping** This product is provided as lyophilized powder which is shipped with ice packs.

Formulation Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH

8.0.


Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as

protectants before lyophilization.

Please refer to the specific buffer information in the printed manual.

**Reconstitution** Please refer to the printed manual for detailed information.

# Data



> 95 % as determined by reducing SDS-PAGE.

## **Background**

Brain-type Natriuretic Peptide (BNP) is a nonglycosylated peptide that is produced predominantly by ventricular

#### For Research Use Only

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com

# **Elabscience Bionovation Inc.**



A Reliable Research Partner in Life Science and Medicine

myocytes and belongs to the natriuretic peptide family. Proteolytic cleavage of the 12 kDa BNP precursor gives rise to Nterminal Pro BNP (NT-proBNP) and mature BNP. N-terminal proB-type natriuretic peptide (NT-proBNP); a useful marker of heart failure (HF); is considered to be secreted mainly from the ventricle; increased serum NT-proBNP levels are also encountered in conditions such as atrial fibrillation (AF) and atrial septal defect in patients without HF.

For Research Use Only

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017 Email: techsupport@elabscience.com

Web: www.elabscience.com