



A Reliable Research Partner in Life Science and Medicine

Biotin Anti-Human CD147 Antibody[HIM6]

Catalog No. E-AB-F1056B Reactivity Human Storage Store at 2~8°C, Avoid freeze / thaw cycles **Applications FCM**

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.

Antigen Information

Alternate Names Basigin,BSG,5F7,Collagenase stimulatory factor,Extracellular matrix metalloproteinase

inducer, EMMPRIN, Leukocyte activation antigen M6,OK blood group antigen, Tumor cell-

derived collagenase stimulatory factor, TCSF, CD147

Uniprot ID P35613

Background CD147, also known as neurothelin or basigin, is a member of the Ig superfamily. It is a 55-65 kD

> type I transmembrane glycoprotein which is primarily expressed on leukocytes, erythrocytes, platelets, and endothelial cells. CD147 is reported to have a function during embryonal brain

development and/or play a role in integrin-mediated adhesion in brain endothelia.

Product Details

Form Liquid Concentration 0.5 mg/mL Size $25 \mu g / 100 \mu g$ Clone No. HIM6 Host Mouse

Mouse IgG1, κ **Isotype**

Reactivity Human **Application FCM**

Isotype Control Biotin Mouse IgG1, κ Isotype Control[MOPC-21] [Product E-AB-F09793B]

Storage Buffer Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.

Biological ice pack at 4 °C Shipping Stability & Storage Keep as concentrated solution. Store at 2~8°C .Do not freeze.

This product is guaranteed up to one year from purchase.

For Research Use Only

Toll-free: 1-888-852-8623 Email: techsupport@elabscience.com Web: www.elabscience.com

Tel: 1-832-243-6086 Fax: 1-832-243-6017

Elabscience Bionovation Inc.



A Reliable Research Partner in Life Science and Medicine

Recommended usage

Each lot of this antibody is quality control tested by flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is $\leq 1.0 \,\mu g$ per 10^6 cells in $100 \,\mu L$ volume or $100 \,\mu L$ of whole blood. It is recommended that the reagent be titrated for optimal performance for each application.

Related Information

- 1. Sample Preparation for Flow Cytometry https://www.elabscience.com/List-detail-5594.html
- 2. Staining Cell Surface Targets for Flow Cytometry https://www.elabscience.com/List-detail-5568.html
- 3. Flow Cytometry Troubleshooting Tips https://www.elabscience.com/List-detail-5593.html
- 4. How to select the appropriate detection channel through the spectrogram? https://www.elabscience.com/List-detail-459742.html

For Research Use Only

Web: www.elabscience.com Email: techsupport@elabscience.com

Toll-free: 1-888-852-8623 Tel: 1-832-243-6086 Fax: 1-832-243-6017