

Recombinant Human AMPK (G1/B1/A1) Heterotrimer Protein



Catalog Number:PKSH030314 1 Publications

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	AMPK;AMPK α 1
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Met 1-Pro 331&Met 1-Ile 270&Met 1-Gln 559
Accession	P54619-1&Q9Y478-1&Q13131-1
Calculated Molecular Weight	38&30&92 kDa
Observed molecular weight	38&40&95 kDa
Tag	N-His-GST
Bioactivity	The specific activity was determined to be 6 nmol/min/mg using synthetic SAMS peptide (HMRSAMSGHLHLVKRR) as substrate.

Properties

Purity	> 94 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μ g of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as sterile solution of 50mM Tris, 200mM NaCl, 2mM GSH, 1mM EDTA, 1mM DTT, 10% glycerol, pH 7.4
Reconstitution	Not Applicable

Data



> 94 % as determined by reducing SDS-PAGE.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017