

Recombinant Human FSTL1 Protein (His Tag, E. coli)

Catalog Number:PKSH032455



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Follistatin-Related Protein 1;Follistatin-Like Protein 1;FSTL1;FRP
Species	Human
Expression Host	E.coli
Sequence	Glu21-Ile308
Accession	Q12841
Calculated Molecular Weight	33.8 kDa
Observed molecular weight	30-38 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Follistatin-Related Protein 1 (FSTL1) is a secreted protein that contains two EF-hand domains; one follistatin-like domain; one Kazal-like domain; and one VWFC domain. Its functional significance in physiological and pathological processes is not completely understood. However, FSTL1 is thought to modulate the action of some growth factors on cell proliferation and differentiation. FSTL1 maybe an autoantigen associated with rheumatoid arthritis.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017