

# Goat Anti-Rabbit IgG

Catalog Number: E-AB-1034

15 Publications



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	Rabbit
<b>Immunogen</b>	Rabbit IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS, PH7.4

## Applications Recommended Dilution

<b>Other</b>	When the antibody is labeled, for example, with HRP, FITC, BIOTIN, etc, BSA and NaN <sub>3</sub> cannot be contained in reagents.
--------------	---

## Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017