

Goat Anti-Rabbit IgG(H+L)(Rhodamine/TRITC conjugated)



Catalog Number:E-AB-1053

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rabbit
Immunogen	Rabbit IgG
Host	Goat
Purification	Affinity purification
Conjugation	TRITC (Ex/Em:550nm/620nm)
Formulation	PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Immunofluorescence	1:20-100
---------------------------	----------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017