

Goat Anti-Mouse IgG(H+L)(peroxidase/HRP conjugated)



Catalog Number: E-AB-1001 55 Publications

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	HRP
Formulation	PBS with 50% glycerol, PH7.4

Applications Recommended Dilution

Western Blotting	1:5000-50000 (Enhanced chemiluminescent detection)
Western Blotting	1:1000-5000(DAB detection)
Direct ELISA	1:5000-30000(TMB detection)
Immunohisto/cytochemistry, 1:400-4000(DAB detection)	

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017