

PPP4R1L Polyclonal Antibody

Catalog No. E-AB-32650

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Monkey
Immunogen	Synthesized peptide derived from the Internal region of human PPP4R1L
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide,0.5% protective protein and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:500-1:2000
ELISA	1:10000

Data



Western Blot analysis of COLO205 cells with
PPP4R1L Polyclonal Antibody.

Observed Mw:50kDa
Calculated Mw:50kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

PPP4R1L (Protein Phosphatase 4 Regulatory Subunit 1 Like (Pseudogene)) is a Pseudogene.

For Research Use Only