

# Recombinant Human RTN4/NOGO-A Protein (GST Tag)



Catalog Number:PKSH030757

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Synonyms</b>	ASY;DGU;DKFZp781L1143;HIGM4;Nbla00271;Nbla10545;NI220/250;NOGO;NOGO-A;Nogo-B;Nogo-C;NOGOC;NSP;NSP-CL;RTN-X;RTN4-A;RTN4-B1;RTN4-B2;RTN4-C;UDG;UNG1;UNG15;UNG2
<b>Species</b>	Human
<b>Expression Host</b>	E.coli
<b>Sequence</b>	Met 1-Val 185
<b>Accession</b>	NP_065393.1
<b>Calculated Molecular Weight</b>	46.2 kDa
<b>Observed molecular weight</b>	48 kDa
<b>Tag</b>	N-GST

## Properties

<b>Purity</b>	> 70 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	Please contact us for more information.
<b>Storage</b>	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
<b>Shipping</b>	This product is provided as lyophilized powder which is shipped with ice packs.
<b>Formulation</b>	Lyophilized from sterile 20mM Tris, 150mM NaCl, 1mM DTT, 0.2mM GSH, pH 7.0 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manu
<b>Reconstitution</b>	Please refer to the printed manual for detailed information.

## Data



> 70 % as determined by reducing SDS-PAGE.

## Background

Reticulon-4, also known as Foccen, Neurite outgrowth inhibitor, Nogo protein, Neuroendocrine-specific protein, Neuroendocrine-specific protein C homolog, RTN-x, Reticulon-5 and RTN4, is a multi-pass membrane protein which contains onereticulon domain. Isoform1of RTN4 is specifically expressed in brain and testis and weakly in heart and skeletal muscle. Isoform2of RTN4 is widely expressed except for the liver. Isoform3of RTN4 is expressed in brain, skeletal muscle and adipocytes. Isoform4of RTN4 is testis-specific. Reticulon-4 / RTN4 is a developmental neurite

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017

# Recombinant Human RTN4/NOGO-A Protein (GST Tag)

Catalog Number:PKSH030757



growth regulatory factor with a role as a negative regulator of axon-axon adhesion and growth, and as a facilitator of neurite branching. Reticulon-4 / RTN4 regulates neurite fasciculation, branching and extension in the developing nervous system. Reticulon-4 / RTN4 is involved in down-regulation of growth, stabilization of wiring and restriction of plasticity in the adult CNS. It regulates the radial migration of cortical neurons via an RTN4R-LINGO1 containing receptor complex. Isoform 2 of RTN4 reduces the anti-apoptotic activity of Bcl-xl and Bcl-2. This is likely consecutive to their change in subcellular location, from the mitochondria to the endoplasmic reticulum, after binding and sequestration. Isoform 2 and isoform 3 of RTN4 inhibit BACE1 activity and amyloid precursor protein processing.

---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: [www.elabscience.com](http://www.elabscience.com)

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)