

Mouse IgG1 Isotype Antibody(Biotin conjugated)

Catalog No.	E-AB-1207	Conjugation	Biotin
Storage	Store at 4°C, Avoid freeze / thaw cycles	Applications	FCM

Important Note:

Centrifuge before opening to ensure complete recovery of vial contents.
This product is guaranteed up to one year from purchase.

Product Information

Conjugation:Biotin

Buffer:PBS with 0.09% NaN₃, 0.2% BSA, PH7.2

Storage:Keep as concentrated solution. Store at 4°C and protected from prolonged exposure to light. Do not freeze.

Usage:Each lot of this antibody is quality control tested by immunofluorescentstaining with flow cytometric analysis (The amount of the reagent issuggested to be used ≤ 0.5 µg/10⁶ cells in 100 µl). Since applicationsvary, the appropriate dilutions must be determined for individual use.

Product Characterization

Clone No.	MOPC-21
Host	Mouse
Isotype	IgG1

Background

The MOPC-21 immunoglobulin has unknown specificity. The isotype of this antibody is mouse IgG1. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

For Research Use Only