

# Anti-Ebola virus EBOV(subtype Sudan, strain Gulu) VP40/Matrix protein VP40 Polyclonal Antibody



Catalog Number:E-AB-V1117

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Reactivity</b>	Ebola virus
<b>Immunogen</b>	A synthetic peptide corresponding to the C-terminus of the Zaire ebolavirus (Guinea 2014) VP40
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Antigen Affinity
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	0.2 µm filtered solution in PBS

## Applications Recommended Dilution

<b>WB</b>	1:1000-1:5000
-----------	---------------

## Data



Western Blot analysis of Recombinant EBOV (subtype Sudan, strain Gulu) VP40 / Matrix protein VP40 Protein (His & MBP Tag)(PKSV030147 with 50ng) using Anti-Ebola virus EBOV(subtype Sudan, strain Gulu) VP40/Matrix protein VP40 Polyclonal Antibody at dilution of 1:1000.

## Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017