

Anti-Ebola virus EBOV(subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 Polyclonal Antibody



Catalog Number: E-AB-V1124

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Ebola virus
Immunogen	Recombinant EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 Protein (His Tag)
Host	Rabbit
Isotype	IgG
Purification	Protein A Affinity
Conjugation	Unconjugated
Formulation	0.2 µm filtered solution in PBS

Applications Recommended Dilution

ELISA	1:1000-1:2000
--------------	---------------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017