

# Bcl-2 Polyclonal Antibody

Catalog Number:E-AB-40526



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Reactivity</b>	Zebrafish
<b>Immunogen</b>	Synthetic peptide corresponding to zebrafish BCL-2 conjugated to keyhole limpet haemocyanin.
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Antigen Affinity Purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS with 0.02% sodium azide, 50% glycerol,pH 7.4

## Applications Recommended Dilution

<b>WB</b>	1:1000-2000
-----------	-------------

## Data



Western blot of Zebrafish whole lysates with anti-BCL-2 rabbit polyclonal antibody at dilution of 1:1000.

**Observed Mw:26kDa**  
**Calculated Mw:26kDa**

## Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

## Background

Bcl-2 (B-cell lymphoma 2), is the founding member of the Bcl-2 family of regulator proteins that regulate cell death (apoptosis), by either inducing (pro-apoptotic) or inhibiting (anti-apoptotic) apoptosis. In mammals, it is found to be associated with

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017