

Recombinant Human CD155/PVR/NECL5 protein (His tag)

Catalog Number:PDMH100147

DIA • AN[®]
by Elabscience

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	Poliovirus Receptor;Nectin-Like Protein 5;NECL-5;CD155;PVR;PVS
Species	Human
Expression Host	HEK293 Cells
Sequence	Met1-Asn343
Accession	P15151
Calculated Molecular Weight	37.6 kDa
Observed molecular weight	60 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Poliovirus Receptor (PVR) is a 70 kDa type I transmembrane single-span glycoprotein that belongs to the nectin-like (Nectin) family and was originally identified based on its ability to mediate the cell attachment and entry of poliovirus (PV); an etiologic agent of the central nervous system disease poliomyelitis. PVR contains three Ig-like extracellular domains; a transmembrane segment; and a cytoplasmic tail. The normal cellular function of PVR maybe the involvement of intercellular adhesion between epithelial cells. Alternate splicing of the PVR mRNA yields four different isoforms (α ; β ; γ ; and δ) with identical extracellular domains.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017