

Goat IgG

Catalog Number:E-AB-1123



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Host	Goat
Isotype	IgG
Purification	Antigen Affinity Purification
Conjugation	Unconjugated
Formulation	0.02M PBS,pH7.2,50% glycerol

Preparation & Storage

Storage	Store at -20°C. Avoid freeze/thaw cycles.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017