

TTR Polyclonal Antibody

Catalog Number:D-AB-10390L



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Rat
Immunogen	Recombinant Monkey TTR protein expressed by Mammalian
Host	Rabbit
Isotype	IgG
Purification	Antigen Affinity Purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide,50% glycerol pH 7.4


Applications Recommended Dilution

WB	1:500-1:1000
IF	1:50-1:200

Data



Western blot with TTR Polyclonal antibody at dilution of 1:1000.lane 1:Human plasma
Observed Mw:16kDa
Calculated Mw:16kDa



Immunofluorescence analysis of C6 cells using TTR Polyclonal Antibody at dilution of 1:200

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017