

Goat Anti-Human IgG(H+L)(Alkaline Phosphatase/AP conjugated)

Catalog No. E-AB-1046

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen: Human IgG

Purification Method: Affinity purification

Conjugation: AP

Concentration: 0.2 mg/mL

Dilution: Western Blotting: 1:1000-3000(Enhanced chemiluminescent detection)

Direct ELISA: 1:1000-3000(TMB detection)

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

Buffer: PBS with 50% glycerol, PH7.4

Storage: Store at -20°C. Avoid freeze / thaw cycles.