


ATXN3 Polyclonal Antibody


Catalog No.	E-AB-40246	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Images



Western Blot analysis of Mouse brain tissue using ATXN3 Polyclonal Antibody at dilution of 1:500



Western Blot analysis of Rat brain tissue using ATXN3 Polyclonal Antibody at dilution of 1:500



Immunohistochemistry of paraffin-embedded Human lung tissue using ATXN3 Polyclonal Antibody at dilution of 1:200

Immunogen Information

Immunogen	Recombinant Human Ataxin-3 protien
Swissprot	P54252
Synonyms	AT3,Ataxin 3,ataxin 3 variant h,ataxin 3 variant m,ataxin 3 variant ref,ATX3,EC 3.4.22,JOS,Machado Joseph disease,MJD, OTTHUMP00000221583,Rsca3

Product Information

Calculated MW	41kDa
Observed MW	41kDa,22kDa
Buffer	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.
Purify	Affinity purification
Dilution	WB 1:500-1:2000 IHC 1:50-1:200

Background

Machado-Joseph disease, also known as spinocerebellar ataxia-3, is an autosomal dominant neurologic disorder. The protein encoded by this gene contains (CAG)_n repeats in the coding region, and the expansion of these repeats from the normal 12-44 to 52-86 is one cause of Machado-Joseph disease. There is a negative correlation between the age of onset and CAG repeat numbers. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.