

Recombinant Human PFK1/PFKM Protein (His & GST Tag)



Catalog Number:PKSH030320

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	6-phosphofructokinase;muscle type;Phosphofructo-1-kinase isozyme A;Phosphofructokinase 1;Phosphohexokinase;PFKM;PFKX;ATP-PFK;GSD7;PFK-1;PFK1;PFKA;PPP1R122
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Thr 2-Val 780
Accession	P08237-1
Calculated Molecular Weight	112.9 kDa
Observed molecular weight	113 kDa
Tag	N-His-GST

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as sterile solution of 20mM Tris, 500mM NaCl, pH 8.5, 10% glycerol
Reconstitution	Not Applicable

Data



> 90 % as determined by reducing SDS-PAGE.

Background

PFK1, also known as PFKM, is a regulatory glycolytic enzyme. PFK1 converts fructose 6-phosphate and ATP into fructose 1,6-bisphosphate (through PFK-1), fructose 2,6-bisphosphate (through PFK-2) and ADP. It is a muscle-type isozyme. There are three phosphofructokinase isozymes in humans: muscle, liver and platelet. These isozymes function as subunits of the mammalian tetramer phosphofructokinase, which catalyzes the phosphorylation of fructose-6-phosphate to fructose-1,6-bisphosphate. Mutations in PFK1 gene have been related with glycogen storage disease type VII, also identified as Tarui disease.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017