

Recombinant Human FES Kinase/Feline sarcoma oncogene Protein (His & GST Tag)



Catalog Number:PKSH030338

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	FPS
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Met 1-Arg 822
Accession	P07332-1
Calculated Molecular Weight	121 kDa
Observed molecular weight	110 kDa
Tag	N-His-GST
Bioactivity	The specific activity was determined to be 200 nmol/min/mg using Poly(Glu:Tyr) 4:1 as substrate.

Properties

Purity	> 75 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as sterile solution of 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol
Reconstitution	Not Applicable

Data



> 75 % as determined by reducing SDS-PAGE.

Background

Proto-oncogene tyrosine-protein kinase Fes/Fps, also known as Proto-oncogene c-Fes, Proto-oncogene c-Fps, Feline sarcoma oncogene, FES and FPS, is a protein which contains oneFCH domain, oneprotein kinase domain and oneSH2 domain. FES is a non-receptor protein tyrosine kinase expressed in hematopoietic progenitors and differentiated myeloid cells. FES is observed in the nuclear, granular and plasma membrane fractions of primary human neutrophils and the myeloid leukemia cell line, HL-60. The nuclear localization is confirmed by immunocytochemistry of neutrophils. FES has been implicated in granulocyte-macrophage colony-stimulating factor (GM-CSF), interleukin-3 (IL-3) and erythropoietin signal transduction. FES has tyrosine-specific protein kinase activity and that activity is required for maintenance of cellular transformation. FES is also involved in normal hematopoiesis. Its chromosomal location has linked it to a specific translocation event identified in patients with acute promyelocytic leukemia.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017