

Recombinant Human NCF2/P67phox Protein (His & GST Tag)



Catalog Number:PKSH030528

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	NCF-2;NOXA2;P67-PHOX;P67PHOX
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Met 1-Val526
Accession	AAH01606.1
Calculated Molecular Weight	87.6 kDa
Observed molecular weight	80-90 kDa
Tag	N-His-GST

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

NCF2, also known as NCF-2 and p67phox, is a subunit of the multi-protein NADPH oxidase complex. NCF2, NCF1, and a membrane bound cytochrome b558 are required for activation of the latent NADPH oxidase. This oxidase produces a burst of superoxide which is delivered to the lumen of the neutrophil phagosome. Mutations in NCF2 gene, as well as in other NADPH oxidase subunits, can result in chronic granulomatous disease, a disease that causes recurrent infections by catalase-positive organisms.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017