

Recombinant Human IL17BR/IL17RB Protein (Fc Tag)

Catalog Number:PKSH030736



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	CRL4;EVI27;IL17BR;IL17RH1
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Gly289
Accession	NP_061195.2
Calculated Molecular Weight	56.9 kDa
Observed molecular weight	67 kDa
Tag	C-hFc
Bioactivity	Immobilized human IL17BR-Fc at 10µg/mL (100µL/well) can bind biotinylated human Fc-IL25. The EC50 of human Fc-IL25 is 4. 5-10.4ng/mL.

Properties

Purity	> 85 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 85 % as determined by reducing SDS-PAGE.

Background

IL17RB (Interleukin 17 Receptor B) is a Protein Coding gene. IL17RB is the receptor for IL17E, the only member of the IL17 family promoting Th2 reactions. IL17RB is induced on human macrophages by IL4 and enhanced by TGFbeta. Overexpression of IL17RB is associated with poor prognosis and the short survival of breast cancer patients. IL17RB/IL17B signaling triggers a substantial increase in cell growth, proliferation, and migration through the activation of NF-kappaB as well as the up-regulation of the Bcl-2. IL17RB may be the only gene expressed in CD4+ T cells whose transcript measurement is correlated with the variation in IgE level in asthmatics. Diseases associated with

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human IL17BR/IL17RB Protein (Fc Tag)

Catalog Number:PKSH030736



IL17RB include Chronic Mucocutaneous Candidiasis and Seborrheic Infantile Dermatitis.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017