

# Recombinant Human NBL1/DAND1 Protein (His Tag)

Catalog Number:PKSH031808



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Synonyms</b>	D1S1733E;DAN;DAND1;NB;N03
<b>Species</b>	Human
<b>Expression Host</b>	HEK293 Cells
<b>Sequence</b>	Met 1-Asp 180
<b>Accession</b>	NP_005371.1
<b>Calculated Molecular Weight</b>	19.3 kDa
<b>Observed molecular weight</b>	27 kDa
<b>Tag</b>	C-His
<b>Bioactivity</b>	Measured by its ability to inhibit BMP4-induced activity in MC3T3-E1 Mouse osteoblastic cells. The ED50 for this effect is typically 5-20 µg/ml in the presence of 50 ng/mL of recombinant human BMP4.

## Properties

<b>Purity</b>	> 98 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Storage</b>	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
<b>Shipping</b>	This product is provided as lyophilized powder which is shipped with ice packs.
<b>Formulation</b>	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
<b>Reconstitution</b>	Please refer to the printed manual for detailed information.

## Data



> 98 % as determined by reducing SDS-PAGE.

## Background

The Dan (Differential screening-selected gene aberrative in neuroblastoma; also known as N03) gene was first identified as the putative rat tumor suppressor gene and encodes a protein structurally related to Cerberus and Gremlin in vertebrates. It is a founding member of the DAN family of secreted proteins; acts as an inhibitor of cell cycle progression and is closely involved in retinoic acid-induced neuroblastoma differentiation. There are at least five mammalian protein members in the evolutionarily conserved Dan family including DAN; Gremlin/DRM; Cer1 (Cerberus-related); Dante and

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017

# Recombinant Human NBL1/DAND1 Protein (His Tag)

Catalog Number:PKSH031808



PRDC (protein related to DAN and cereberus); and share the C-terminal cystine-knot motif. As a secreted glycoprotein; DAN is a member of a class of glycoproteins shown to be secreted inhibitors of the transforming growth factor-beta (TGF-beta) and bone morphogenic protein pathways. It binds to BMPs and preventing their interactions with signaling receptor complexes; and accordingly regulates the processes of embryonic development and tissue differentiation. DAN gene product may have an important role in regulation of the entry of cells into the S phase. In addition; DAN gene product possesses an ability to revert phenotypes of transformed rat fibroblasts and represents a candidate tumour suppressor gene for neuroblastoma.

---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: [www.elabscience.com](http://www.elabscience.com)

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)