

Recombinant Human NCAM1 Protein (ECD, His Tag)

Catalog No. PKSH031994

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	CD56;NCAM-1;CD56 antigen;MSK39;N-CAM-1;NCAM-1;neural cell adhesion molecule 1;neural cell adhesion molecule;NCAM
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Pro603
Accession	NP_001070150.1
Calculated Molecular Weight	66.1 kDa
Tag	C-His
Bioactivity	Not validated for activity

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

NCAM1, also known as CD56, is a neural adhesion protein (NCAM) which belongs to the immunoglobulin superfamily. NCAM is involved in neural development and in plasticity in the adult brain. UCHL1 is a novel interaction partner of both

For Research Use Only

NCAM isoforms that regulates their ubiquitination and intracellular trafficking. NCAM1 is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc. NCAM1 has also been shown to be involved in the expansion of T cells and dendritic cells which play an important role in immune surveillance.