

Recombinant Human AGR3 Protein (His Tag)

Catalog Number:PKSH032080



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Anterior Gradient Protein 3 Homolog;AG-3;AG3;hAG-3;Breast Cancer Membrane Protein 11;AGR3;BCMP11
Species	Human
Expression Host	HEK293 Cells
Sequence	Ile22-Leu166
Accession	Q8TD06
Calculated Molecular Weight	18.0 kDa
Observed molecular weight	16 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Glycine-HCl, 10% Trehalose, 0.05% Tween80, pH3.5. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer info
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Anterior Gradient Protein 2(AG-2) and Anterior Gradient Protein 3 (AG-3) are human homologues of genes involved in differentiation, are associated with oestrogen receptor-positive breast tumours and interact with metastasis gene C4.4a and dystroglycan (hAG-3 protein). AG-3 could serve as a prognostic marker for survival in patients with low grade and high grade serous ovarian carcinomas.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017