

Recombinant Human Arginase-1/ARG1 Protein (E.coli, His Tag)



Catalog Number:PKSH032091

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Arginase-1;Liver-type arginase;Type I arginase;ARG1
Species	Human
Expression Host	E.coli
Sequence	Met 1-lys322
Accession	P05089
Calculated Molecular Weight	35.8 kDa
Observed molecular weight	40 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 20% Glycerol, 1mM DTT, pH 7.4.
Reconstitution	Not Applicable

Data



> 95 % as determined by reducing SDS-PAGE.

Background

ARG1 is a member of the ureohydrolase family of enzymes. ARG1 can catalyze the hydrolysis of arginine to ornithine and urea. In the urea cycle, ARG1 catalyzes the fifth and final step, a series of biochemical reactions in mammals during which the body disposes of harmful ammonia. ARG1 is a cytosolic enzyme and expressed widely in the liver as part of the urea cycle. Inherited deficiency of this ARG1 causes argininemia, which is an autosomal recessive disorder characterized by hyperammonemia.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017