

Recombinant Human Survivin/BIRC5 Protein

Catalog Number:PKSH032109



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Baculoviral IAP Repeat-Containing Protein 5;Apoptosis Inhibitor 4;Apoptosis Inhibitor Survivin;BIRC5;API4;IAP4
Species	Human
Expression Host	E.coli
Sequence	Met 1-Asp142
Accession	O15392
Calculated Molecular Weight	16.4 kDa
Observed molecular weight	16-18 kDa
Tag	None

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 7.5. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the print
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Baculoviral IAP Repeat-Containing Protein 5 (BIRC5) belongs to the IAP family. BIRC5 exists as a homodimer in the apo state and as a monomer in the CPC-bound state. BIRC5 contains one BIR repeat and is expressed only in fetal kidney and liver, and to lesser extent, lung and brain. BIRC5 functions as a multitasking protein that has dual roles in promoting cell proliferation and preventing apoptosis. BIRC5 is also a component of a chromosome passage protein complex (CPC), which is essential for chromosome alignment and segregation during mitosis and cytokinesis. BIRC5 acts as an important regulator of the localization of this complex.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017