

Recombinant Human CA11 Protein (His Tag)

Catalog Number:PKSH032166



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Carbonic Anhydrase-Related Protein 11;CA-RP XI;CA-XI;CARP XI;Carbonic Anhydrase-Related Protein 2;CA-RP II;CARP-2;CA11;CARP2
Species	Human
Expression Host	HEK293 Cells
Sequence	His24-Arg328
Accession	O75493
Calculated Molecular Weight	35.1 kDa
Observed molecular weight	38 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μ g of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Reconstitution	Not Applicable

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Carbonic Anhydrase-Related Protein 11 (CA11) is a secreted protein member of the α -carbonic anhydrase family. Carbonic Anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. CA11 is expressed abundantly in the brain with moderate expression also present in spinal cord and thyroid. CA11 may play a general role in the central nervous system.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017