

Recombinant Human HSPB2/MKBP Protein (His Tag)

Catalog Number:PKSH032523



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Heat shock protein beta-2;HspB2;DMPK-binding protein;MKBP;
Species	Human
Expression Host	E.coli
Sequence	Met 1-Pro182
Accession	Q16082
Calculated Molecular Weight	21.3 kDa
Observed molecular weight	22-25 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM Tris-HCl, 150mM NaCl, 1mM EDTA, pH 8.0. Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Heat shock protein beta-2(HSPB2) is a protein that in humans is encoded by the HSPB2 gene. HSPB2 belongs to the superfamily of small heat-shock proteins containing a conservative alpha-crystallin domain at the C-terminal part of the molecule. It is expressed preferentially in the heart and skeletal muscle. HSPB2 has been shown to interact with TRAF6, HSPB8, Myotonic dystrophy protein kinase and CRYAB. HSPB2 regulates Myotonic Dystrophy Protein Kinase, which plays an important role in maintenance of muscle structure and function.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017