

Recombinant Human HMGB1 Protein

Catalog Number:PKSH032544



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	High Mobility Group Protein B1;High Mobility Group Protein 1;HMG-1;HMGB1;HMG1
Species	Human
Expression Host	E.coli
Sequence	Pro92-Val176
Accession	P09429
Calculated Molecular Weight	10.0 kDa
Observed molecular weight	13 kDa
Tag	N-MARI

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM HEPES-Na, 500mM NaCl, 0.6mM DTT, pH 7.9. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information i
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

High mobility group protein B1 is a member of the HMGB family consisting of three members; HMGB1; HMGB2 and HMGB3. It contains 2 HMG box DNA-binding domains entitled box A and box B and it is a highly negative-charged C terminus. As a nuclear protein; HMGB1 stabilizes nucleosomes and allows bending of DNA that facilitates gene transcription which is essential for individual survival. Meanwhile; it is revealed that HMGB1 can also act as a cytokine extracellularly and regulates monocyte; T cell; dendritic cell activities in inflammatory responses.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017