

Recombinant Human HMGB3 Protein (His Tag)

Catalog Number:PKSH032546



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	High Mobility Group Protein B3;High Mobility Group Protein 2a;HMG-2a;High Mobility Group Protein 4;HMG-4;HMGB3;HMG2A;HMG4
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Glu200
Accession	O15347
Calculated Molecular Weight	24.0 kDa
Observed molecular weight	26 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

High Mobility Group Protein B3 (HMGB3) belongs to the HMGB family. Members of the HMG box subfamily are thought to have an important role in DNA replication, nucleosome assembly and transcription. HMGB3 binds preferentially single-stranded DNA and unwinds double stranded DNA. HMGB3 consists of 200 amino acids and is localized to the cell nucleus. It contains two HMG box DNA-binding domain. HMGB3 binds preferentially single-stranded DNA and unwinds double stranded DNA.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017