

Recombinant Human MYDGF Protein (His Tag)

Catalog No. PKSH032635

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	UPF0556 protein C19orf10;stromal cell-derived growth factor SF20;C19orf10;Myeloid-derived growth factor;MYDGF
Species	Human
Expression Host	E.coli
Sequence	Ser33-Leu173
Accession	Q969H8
Calculated Molecular Weight	18.0 kDa
Observed molecular weight	17 kDa
Tag	N-His
Bioactivity	Not validated for activity

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 4mM HCl. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Myeloid-derived growth factor (MYDGF) is a secreted protein which belongs to the UPF0556 family. MYDGF was

For Research Use Only

strongly expressed in spleen, prostate and lung, and weakly expressed in the left ventricle and liver. Bone marrow-derived monocyte and paracrine-acting protein promotes cardiac myocyte survival and adaptive angiogenesis for cardiac protection and/or repair after myocardial infarction (MI). MYDGF stimulates endothelial cell proliferation through a MAPK1/3-, STAT3- and CCND1-mediated signaling pathway. It inhibits cardiac myocyte apoptosis in a PI3K/AKT-dependent signaling pathway. MYDGF is involved in endothelial cell proliferation and angiogenesis. It may serve as a prototypical example for the development of protein-based therapies for ischemic tissue repair.