

Recombinant Human NLGN4X Protein (His Tag)

Catalog Number:PKSH032797



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Neurologin-4 X-Linked;Neurologin X;HNLX;NLGN4X;KIAA1260;NLGN4
Species	Human
Expression Host	HEK293 Cells
Sequence	Gln42-Ser676
Accession	Q8N0W4
Calculated Molecular Weight	72.5 kDa
Observed molecular weight	95 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Neurologin 4, X-Linked (NLGN4X) is a single-pass type I membrane protein that belongs to the type-B carboxylesterase/lipase family. NLGN4X is detected at higher levels in heart and at lower levels in the liver, skeletal muscle, and pancreas. NLGN4X is a putative neuronal cell surface protein involved in cell-cell-interactions. NLGN4X may act as splice site-specific ligands for β -neurexins. It has been shown that NLGN4X is involved in the formation and remodeling of central nervous system synapses. NLGN4X also interacts with discs, large (Drosophila) homolog 4 (DLG4). Defects in NLGN4X have been associated with autism and Asperger syndrome.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017