

Recombinant Human NDRG1 Protein (His Tag)

Catalog Number:PKSH032818



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Protein NDRG1;Differentiation-Related Gene 1 Protein;DRG-1;N-myc Downstream-Regulated Gene 1 Protein;Nickel-Specific Induction Protein Cap43;Reducing Agents and Tunicamycin-Responsive Protein;RTP;Rit42;NDRG1;CAP43;DRG1;RTP
Species	Human
Expression Host	E.coli
Sequence	Met 1-Cys394
Accession	Q92597
Calculated Molecular Weight	45.0 kDa
Observed molecular weight	52 kDa
Tag	N-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

Protein NDRG1 is a member of the N-Myc Downregulated Gene family, which is part of the α/β Hydrolase superfamily. Protein NDRG1 is a cytoplasmic protein that is involved in stress responses, hormone responses, cell growth and differentiation. Protein NDRG1 is necessary for p53-mediated caspase activation and apoptosis. Protein NDRG1 mutations are reported to be the cause for hereditary motor and sensory neuropathy-Lom. Decreased NDRG1 expression in glioma is linked to tumor progression; overexpression of NDRG1 is connected to malignant status of esophageal

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human NDRG1 Protein (His Tag)

Catalog Number:PKSH032818



cancer.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com