

Recombinant Human Sclerostin/SOST protein (His tag)

Catalog No. PKSH033010

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	CDD, DAND6, SOST1, VBCH
Species	Human
Expression Host	HEK293 Cells
Sequence	Met1-Tyr213
Accession	Q9BQB4
Calculated Molecular Weight	22.3 kDa
Observed molecular weight	35 kDa
Tag	C-His
Bioactivity	Testing in progress

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Sclerostin, also known as SOST, is a member of the Cerberus/DAN family of BMP antagonists. SOST is a secreted glycoprotein with a C-terminal cysteine knot-like (CTCK) domain. It shows sequence similarity to the DAN (differential

For Research Use Only

screening-selected gene aberrative in neuroblastoma) family of bone morphogenetic protein (BMP) antagonists. SOST is produced by the osteocyte and has anti-anabolic effects on bone formation. It is a negative regulator of bone growth that acts through inhibition of Wnt signaling and bone formation.