

Recombinant Human SOAT2/ACAT2 Protein (His&Trx Tag)



Catalog Number:PKSH033082

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Acetyl-CoA acetyltransferase cytosolic;Acetyl-CoA transferase-like protein;Cytosolic acetoacetyl-CoA thiolase;ACTL;ACAT2
Species	Human
Expression Host	E.coli
Sequence	Met 1-Glu397
Accession	Q9BWD1
Calculated Molecular Weight	59.1 kDa
Observed molecular weight	60 kDa
Tag	N-Trx-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
Reconstitution	Not Applicable

Data



> 95 % as determined by reducing SDS-PAGE.

Background

ACAT2 is a cytoplasmic enzyme which belongs to the thiolase family. ACAT2 takes part in lipid metabolism, lipoprotein assembly, catalyzing cholesterol esterification in mammalian cells. It is responsible for the synthesis of cholesteryl esters which are part of lipoproteins containing apoB. ACAT2 deficiency contributes to severe mental retardation and hypotonus.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017