

Recombinant Human Syntaxin-6/STX6 Protein

Catalog Number:PKSH033092



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	Syntaxin-6;STX6
Species	Human
Expression Host	E.coli
Sequence	Ser2-Gln234
Accession	O43752
Calculated Molecular Weight	26.8 kDa
Observed molecular weight	32 kDa
Tag	None

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual
Reconstitution	Please refer to the printed manual for detailed information.

Background

Syntaxin-6 (STX6) is a single-pass type IV membrane protein, which belongs to the Syntaxin family. STX6 is mainly localized in the plasma membrane. STX6 contains one t-SNARE coiled-coil homology domain and involved in intracellular vesicle trafficking. When STX6 function is inhibited, internalization through caveolae is dramatically reduced, whereas other endocytic mechanisms are unaffected. It is reported that STX6 is necessary for proper expression of focal adhesion kinase and integrin alpha5 adhesion receptor.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017