

Recombinant Human Thioredoxin-2/TXN2 Protein

Catalog Number:PKSH033109



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Thioredoxin Mitochondrial;MTRX;Mt-Trx;Thioredoxin-2;TXN2;TRX2
Species	Human
Expression Host	E.coli
Sequence	Thr60-Gly166
Accession	Q99757
Calculated Molecular Weight	12.0 kDa
Observed molecular weight	13 kDa
Tag	None

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Thioredoxin-2 (TXN2) is a mitochondrial member of the thioredoxin family. Thioredoxin-2 is extensively expressed in adult and fetal tissues. Thioredoxin-2 contains an N-terminal 59 amino acid transit peptide; which is cleaved before translocating to mitochondria. Mitochondrial thioredoxin play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis. Thioredoxin-2 could be involved in the resistance to anti-tumor agents and possesses a dithiol-reducing activity. In addition; Thioredoxin-2 is important at low oxidative stress conditions.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017