

Recombinant Human VMO1 Protein (His Tag)

Catalog Number:PKSH033217



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Vitelline Membrane Outer Layer Protein 1 Homolog;VMO1
Species	Human
Expression Host	HEK293 Cells
Sequence	Gln25-Ser202
Accession	Q7Z5L0
Calculated Molecular Weight	20.1 kDa
Observed molecular weight	18-22 kDa
Tag	C-His

Properties

Purity	> 85 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, 0.5mM EDTA, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 85 % as determined by reducing SDS-PAGE.

Background

Vitelline membrane outer layer protein 1 homolog (VMO1) belongs to the VMO1 family is a 202 amino acid secreted protein. Exact function not known, component of the outer membrane of the vitelline layer of the egg. Seems to be able to synthesize N-acetylchito-oligosaccharides (n=14-15) from hexasaccharides of N-acetylglucosamine in a manner similar to the transferase activity of lysozyme.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017