

Recombinant Mouse Galectin 9/LGALS9 Protein (GST Tag)

Catalog No. PKSM041336

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Galectin-9;Lgals9;Gal-9
Species	Mouse
Expression Host	E.coli
Sequence	Met1-Thr322
Accession	O08573-2
Calculated Molecular Weight	62.8 kDa
Observed molecular weight	58 kDa
Tag	N-GST
Bioactivity	Not validated for activity

Properties

Purity	> 85 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 1mM EDTA, 5mM β-Mercaptoethanol, 5% Trehalose, 400mM NaCl, pH6.5.
Reconstitution	Not Applicable

Data



> 85 % as determined by reducing SDS-PAGE.

Background

Galectin-9 is a cytoplasmic protein that contains two galectin domains. Galectin-9 is an S-type lectin that is over-expressed in Hodgkin's disease tissue. Galectin-9 binds galactosides and has high affinity for the Forssman pentasaccharide. Galectin-9 plays a role in thymocyte-epithelial interactions relevant to the biology of the thymus and inhibits cell proliferation. Galectin-9 is a ligand for HAVCR2/TIM3 and induces T-helper type 1 lymphocyte (Th1) death. In addition, Galectin-9 suppresses tumor cell metastasis by interfering with the associations CD44, VCAM-1, Integrin

For Research Use Only

$\alpha 4\beta 1$.

For Research Use Only

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017