

# Recombinant Vaccinia virus Soluble Interferon $\alpha/\beta$ Receptor B18 Protein (His Tag)

Catalog Number:PKSQ050089

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Synonyms</b>	Soluble interferon alpha/beta receptor B18;VACWR200
<b>Species</b>	Vaccinia virus
<b>Expression Host</b>	HEK293 Cells
<b>Sequence</b>	His20-Glu351
<b>Accession</b>	P25213
<b>Calculated Molecular Weight</b>	39.2 kDa
<b>Observed molecular weight</b>	50-65 kDa
<b>Tag</b>	C-His

## Properties

<b>Purity</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	< 1.0 EU per $\mu\text{g}$ of the protein as determined by the LAL method.
<b>Storage</b>	Generally, lyophilized proteins are stable for up to 12 months when stored at $-20$ to $-80^{\circ}\text{C}$ . Reconstituted protein solution can be stored at $4-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.
<b>Shipping</b>	This product is provided as lyophilized powder which is shipped with ice packs.
<b>Formulation</b>	Lyophilized from a $0.2 \mu\text{m}$ filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
<b>Reconstitution</b>	Please refer to the printed manual for detailed information.

## Data



> 95 % as determined by reducing SDS-PAGE.

## Background

B18R, also known as B19R in the Copenhagen strain of Vaccinia, is a type I interferon (IFN)-binding protein encoded by the B18R open reading frame in the Western Reserve strain of vaccinia virus. B18R exists in two forms, soluble and membrane-bound form and also has a broad species specificity. B18R has high affinity for human IFN-alpha and also binds rabbit, bovine, rat, pig, and mouse IFN-alpha and IFN-beta. Secreted B18R binds to uninfected and infected cells. B18R presents at the cell surface and protects cells from the antiviral state. It has shown that binding of soluble recombinant B18R protects cultured cells from IFN and allows vaccinia virus replication.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017