

REG4 Polyclonal Antibody

Catalog No. E-AB-40314

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant Mouse Regenerating islet-derived protein 4 protein
Host	Rabbit
Isotype	IgG
Purification	Antigen Affinity Purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.


Applications Recommended Dilution

IHC 1:100-1:400


Data




Immunohistochemistry of paraffin-embedded Human colon using REG4 Polyclonal Antibody at dilution of 1:200.



Immunohistochemistry of paraffin-embedded Human duodenum using REG4 Polyclonal Antibody at dilution of 1:200.



Immunohistochemistry of paraffin-embedded Mouse colon using REG4 Polyclonal Antibody at dilution of 1:200 [100x].



Immunohistochemistry of paraffin-embedded Mouse colon using REG4 Polyclonal Antibody at dilution of 1:200 [400x].

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Calcium-independent lectin displaying mannose-binding specificity and able to maintain carbohydrate recognition activity

For Research Use Only

in an acidic environment. May be involved in inflammatory and metaplastic responses of the gastrointestinal epithelium.