

S100A9 Polyclonal Antibody

Catalog No. E-AB-40317

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant Rat Protein S100-A9 protein
Host	Rabbit
Isotype	IgG
Purification	Antigen Affinity Purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide,1% protective protein and 50% glycerol,pH7.4

Applications Recommended Dilution


WB	1:500-1:1000
IHC	1:100-1:300

Data




Western Blot analysis of Rat spleen using S100A9 Polyclonal Antibody at dilution of 1:1000.


Observed Mw:13 kDa
Calculated Mw:13 kDa



Immunohistochemistry of paraffin-embedded Human breast cancer using S100A9 Polyclonal Antibody at dilution of 1:200.




Immunohistochemistry of paraffin-embedded Human kidney using S100A9 Polyclonal Antibody at dilution of 1:200.




Immunohistochemistry of paraffin-embedded Human liver cancer using S100A9 Polyclonal Antibody at dilution of 1:200.

For Research Use Only



Immunohistochemistry of paraffin-embedded Rat lung using S100A9 Polyclonal Antibody at dilution of 1:200.



Immunohistochemistry of paraffin-embedded Rat stomach using S100A9 Polyclonal Antibody at dilution of 1:200. 100x.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21.