

## Goat Anti-Mouse IgG F(ab')<sub>2</sub>(peroxidase/HRP conjugated)

Catalog No. E-AB-1101

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Immunogen</b>	Mouse IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	HRP
<b>Buffer</b>	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

### Applications Recommended Dilution

<b>WB</b>	1:8000-15000 (Enhanced chemiluminescent detection)
<b>WB</b>	1:1000-5000(DAB detection)
<b>Direct ELISA</b>	1:4000-8000(TMB detection)
<b>IHC</b>	1:400-4000(DAB detection)

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.