

Goat Anti-Rabbit IgG F(ab')₂(Elab Fluor® 488 conjugated)

Catalog No. E-AB-1103

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Immunogen	Rabbit IgG
Host	Goat
Purification	Affinity purification
Conjugation	Elab Fluor®488 (Ex/Em:493nm/519nm)
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

Immunofluorescence	1:100-300
---------------------------	-----------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---

For Research Use Only