

## Goat Anti-Mouse IgG F(ab')<sub>2</sub>(Elab Fluor® 488 conjugated)

Catalog No. E-AB-1104

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

### Description

<b>Immunogen</b>	Mouse IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Elab Fluor®488 (Ex/Em:493nm/519nm)
<b>Buffer</b>	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

### Applications Recommended Dilution

<b>Immunofluorescence</b>	1:100-300
---------------------------	-----------

### Preparation & Storage

<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---