

Goat Anti-Rabbit IgG F(ab')₂(peroxidase/HRP conjugated)

Catalog No. E-AB-1102

Note: Centrifuge before opening to ensure complete recovery of vial contents.**Description**

Immunogen	Rabbit IgG
Host	Goat
Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications **Recommended Dilution**

WB	1:8000-15000 (Enhanced chemiluminescent detection)
WB	1:1000-5000(DAB detection)
Direct ELISA	1:4000-8000(TMB detection)
IHC	1:400-4000(DAB detection)

Preparation & Storage**Storage** Store at -20°C. Avoid freeze / thaw cycles.**For Research Use Only**