

# Goat Anti-Rabbit IgG F(ab')<sub>2</sub>(peroxidase/HRP conjugated)

Catalog Number: E-AB-1102

2 Publications

**DIA • AN**<sup>®</sup>  
by Elabscience

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Immunogen</b>	Rabbit IgG
<b>Host</b>	Goat
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	HRP
<b>Formulation</b>	PBS with 50% glycerol, PH7.4

## Applications Recommended Dilution

<b>Western Blotting</b>	1:8000-15000 (Enhanced chemiluminescent detection)
<b>Western Blotting</b>	1:1000-5000(DAB detection)
<b>Direct ELISA</b>	1:4000-8000(TMB detection)
<b>Immunohisto/cytoche mistry</b>	1:400-4000(DAB detection)

## Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017