

Goat Anti-Rabbit IgG(H+L)(peroxidase/HRP conjugated)

Catalog No. E-AB-1003

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Rabbit
Immunogen	Rabbit IgG
Host	Goat
Purification	Affinity purification
Conjugation	HRP
Buffer	PBS with 0.1% PC300, 1% protective protein and 50% glycerol, pH 7.4

Applications Recommended Dilution

WB	1:5000-50000(Enhanced chemiluminescent detection)
WB	1:1000-5000(DAB detection)
Direct ELISA	1:5000-30000(TMB detection)
IHC	1:500-2000(DAB detection)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.