

# Rabbit Anti-Rat IgG(H+L)(FITC conjugated)

Catalog Number:E-AB-1018

1 Publications



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

|                     |                              |
|---------------------|------------------------------|
| <b>Reactivity</b>   | Rat                          |
| <b>Immunogen</b>    | Rat IgG                      |
| <b>Host</b>         | Rabbit                       |
| <b>Purification</b> | Affinity purification        |
| <b>Conjugation</b>  | FITC (Ex/Em:495nm/525nm)     |
| <b>Formulation</b>  | PBS with 50% glycerol, PH7.4 |

## Applications Recommended Dilution

|                           |          |
|---------------------------|----------|
| <b>Immunofluorescence</b> | 1:20-100 |
|---------------------------|----------|

## Preparation & Storage

|                |   |
|----------------|---|
| <b>Storage</b> | Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light. |
|----------------|---|

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017